

# C21orf45 antibody



Catalog Number: 108705

## Product name

C21orf45 antibody

## Immunogen

[Human C21orf45 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

C21orf45 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 105 kDa

## Preparation

This antibody was obtained by immunization of C21orf45 recombinant protein (Accession Number: NM\_018944). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108705(C21orf45 Antibody) at dilution of 1:300



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:108705(C21orf45 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)