

C21orf58 antibody



Catalog Number: 108706

Product name

C21orf58 antibody

Specificity

Human; other species not tested.

Antibody description

C21orf58 Rabbit Polyclonal antibody. Positive IHC detected in human skin tissue, human brain tissue.

Preparation

This antibody was obtained by immunization of C21orf58 recombinant protein (Accession Number: BC121008). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:108706(C21orf58 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:108706(C21orf58 Antibody) at dilution of 1:50 (under 40x lens)