

# C2orf49 antibody

Catalog Number: 108714

---

**Product name**

C2orf49 antibody

**Specificity**

human, mouse, rat; other species not tested.

**Antibody description**

C2orf49 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse liver tissue. Observed molecular weight by Western-blot: 35kd

**Preparation**

This antibody was obtained by immunization of C2orf49 recombinant protein (Accession Number: XM\_047445805). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:200-1:2000

**Validations**

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108714(C2orf49 antibody) at dilution of 1:500