

C2orf60 antibody



Catalog Number: 108716

Product name

C2orf60 antibody

Specificity

Human; other species not tested.

Antibody description

C2orf60 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of C2orf60 recombinant protein (Accession Number: NM_001039693). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108716(C2orf60 Antibody) at dilution of 1:600