

C2orf60 antibody



Catalog Number: 108716

Product name

Polyclonal

C2orf60 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

C2orf60 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 40-45 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of C2orf60 recombinant protein (Accession Number: NM_001039693). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108716(C2orf60 Antibody) at dilution of 1:600

Store at -20°C. DO NOT ALIQUOT

Clonality