

C3orf10 antibody



Catalog Number: 108719

Product name

C3orf10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

C3orf10 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue.

Preparation

This antibody was obtained by immunization of C3orf10 recombinant protein (Accession Number: NM_018462). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:108719(C3orf10 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:108719(C3orf10 antibody) at dilution of 1:50 (under 40x lens)