

C4orf14 antibody



Catalog Number: 108730

Product name

C4orf14 antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

C4orf14 Rabbit Polyclonal antibody. Positive IHC detected in human colon. Positive WB detected in mouse kidney tissue, mouse brain tissue. Observed molecular weight by Western-blot: 85-95 kDa

Preparation

This antibody was obtained by immunization of C4orf14 recombinant protein (Accession Number: NM_032313). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108730(C4orf14 antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human colon using Catalog No:108730(C4orf14 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon using Catalog No:108730(C4orf14 antibody) at dilution of 1:100 (under 40x lens)