

C4orf29 antibody



Catalog Number: 108731

Product name

C4orf29 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

C4orf29 Rabbit Polyclonal antibody. Positive IHC detected in human colon tissue, human tonsillitis tissue. Positive IF detected in A549 cells. Positive WB detected in PC-3 cells, A549 cells, mouse heart tissue. Observed molecular weight by Western-blot: 60 kDa




Preparation

This antibody was obtained by immunization of C4orf29 recombinant protein (Accession Number: NM_001366048). Purification method: Antigen affinity purified.

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108731(C4orf29 Antibody) at dilution of 1:600

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:108731(C4orf29 Antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

Immunohistochemistry of paraffin-embedded

C4orf29 antibody

Catalog Number: 108731

human colon tissue slide using Catalog No:108731(C4orf29 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of A549 cells using Catalog No:108731(C4orf29 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)