

CA1 antibody

Catalog Number: 108741

Product name

CA1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CA1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse spleen tissue, rat spleen tissue. Positive IP detected in K-562 cells. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of CA1 recombinant protein (Accession Number: NM_001164830). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:50-1:500

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108741(CA1 antibody) at dilution of 1:500



IP Result of anti-CA1 (IP:Catalog No:108741, 4ug; Detection:Catalog No:108741 1:500) with K-562 cells lysate 1200ug.