

CA10 antibody

Catalog Number: 108742

Product name

CA10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CA10 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of CA10 recombinant protein (Accession Number: NM_001082534). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:108742(CA10 antibody) at dilution of 1:100 (under 10x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108742(CA10 Antibody) at dilution of 1:700