

# CA10 antibody



Catalog Number: 108742

## Product name

CA10 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CA10 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of CA10 recombinant protein (Accession Number: NM\_001082534). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:108742(CA10 antibody) at dilution of 1:100 (under 10x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108742(CA10 Antibody) at dilution of 1:700