

CA5B antibody



Catalog Number: 108750

Product name

CA5B antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

CA5B Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 32kd

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of CA5B recombinant protein (Accession Number: NM_007220). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108750(CA5B antibody) at dilution of 1:500