

CABC1 antibody



Catalog Number: 108754

Product name

CABC1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CABC1 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human heart tissue, mouse heart tissue. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 60kd

Preparation

This antibody was obtained by immunization of CABC1 recombinant protein (Accession Number: XM_011544240). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:108754(CABC1 antibody) at dilution of 1:800



Immunofluorescent analysis of Hela cells, using CABC1 antibody Catalog No:108754 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).