

# CACNA1G antibody



Catalog Number: 108759

## Product name

CACNA1G antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CACNA1G Rabbit Polyclonal antibody. Positive IF detected in mouse heart tissue.

## Preparation

This antibody was obtained by immunization of CACNA1G recombinant protein (Accession Number: NM\_198382). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF

## Dilutions


Recommended Dilution:

IF: 1:20-1:200

## Validations

CACNA1G Polyclonal Antibody

Immunofluorescence of frozen mouse heart section. 1:100 dilution overnight at 4°C in Block. Secondary antibody: Alexa 594 Goat anti-rabbit. Apparent sarcolemmal expression, as expected.



IF result of anti-CACNA1G (Catalog No: 108759) in frozen mouse heart section by Dr Nicola Smart.