

# CACNB4 antibody



Catalog Number: 108764

## Product name

Polyclonal

CACNB4 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

CACNB4 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Dilutions

## Preparation

Recommended Dilution:


This antibody was obtained by immunization of CACNB4 recombinant protein (Accession Number: NM\_000726). Purification method: Antigen affinity purified.

WB: 1:500-1:5000

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Validations



## Storage

Store at -20°C. DO NOT ALIQUOT

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108764(CACNB4 Antibody) at dilution of 1:1000

## Clonality