

CACNG2 antibody



Catalog Number: 108766

Product name

CACNG2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CACNG2 Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006078).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:108766(CACNG2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:108766(CACNG2 Antibody) at dilution of 1:50 (under 40x lens)