

CACNG5 antibody



Catalog Number: 108768

Product name

Polyclonal

CACNG5 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

CACNG5 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, human brain tissue. Observed molecular weight by Western-blot: 35kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of CACNG5 recombinant protein (Accession Number: NM_145811). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:108768(CACNG5 antibody) at dilution of 1:300