

C7orf49 antibody

Catalog Number: 108790

Product name

C7orf49 antibody

Specificity

Human; other species not tested.

Antibody description

C7orf49 Rabbit Polyclonal antibody. Positive IHC detected in human skin tissue.

Preparation

This antibody was obtained by immunization of C7orf49 recombinant protein (Accession Number: BC000168). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:10-1:100

Validations



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:108790(C7orf49 Antibody) at dilution of 1:25 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:108790(C7orf49 Antibody) at dilution of 1:25 (under 40x lens)