

C9orf86 antibody



Catalog Number: 108805

Product name

C9orf86 antibody

Specificity

Human; other species not tested.

Antibody description

C9orf86 Rabbit Polyclonal antibody. Positive WB detected in MDA-MB-453s cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 125 kDa

Preparation

This antibody was obtained by immunization of C9orf86 recombinant protein (Accession Number: BC002945). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IF: 1:50-1:500

Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:108805(C9orf86 Antibody) at dilution of 1:600



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:108805(C9orf86 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)