

Cadherin-20 antibody



Catalog Number: 108809

Product name

Cadherin-20 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Cadherin-20 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human placenta tissue, mouse brain tissue. Observed molecular weight by Western-blot: 120-130kd

Preparation

This antibody was obtained by immunization of Cadherin-20 recombinant protein (Accession Number: NM_031891). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108809(CDH20 antibody) at dilution of 1:1000