

CADM4 antibody



Catalog Number: 108814

Product name

CADM4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CADM4 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, COLO 320 cells, mouse lung tissue. Observed molecular weight by Western-blot: 65-70 kDa

Preparation

This antibody was obtained by immunization of CADM4 recombinant protein (Accession Number: NM_145296). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108814(CADM4 antibody) at dilution of 1:1000