

Calcyphosine antibody



Catalog Number: 108819

Product name

Calcyphosine antibody

Specificity

Human; other species not tested.

Antibody description

Calcyphosine Rabbit Polyclonal antibody. Positive WB detected in human colon tissue, human spleen tissue. Observed molecular weight by Western-blot: 24 kDa

Preparation

This antibody was obtained by immunization of Calcyphosine recombinant protein (Accession Number: NM_004058). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:108819(CAPS antibody) at dilution of 1:500



human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:108819(CAPS antibody) at dilution of 1:500