

# Calcyphosine antibody



Catalog Number: 108819

## Product name

Calcyphosine antibody

## Specificity

Human; other species not tested.

## Antibody description

Calcyphosine Rabbit Polyclonal antibody. Positive WB detected in human colon tissue, human spleen tissue. Observed molecular weight by Western-blot: 24 kDa

## Preparation

This antibody was obtained by immunization of Calcyphosine recombinant protein (Accession Number: NM\_004058). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:108819(CAPS antibody) at dilution of 1:500



human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:108819(CAPS antibody) at dilution of 1:500