

Product name

CALML3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CALML3 Rabbit Polyclonal antibody. Positive WB detected in A431 cells. Positive IF detected in A431 cells. Observed molecular weight by Western-blot: 17 kDa

Preparation

This antibody was obtained by immunization of CALML3 recombinant protein (Accession Number: NM_005185). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

IF: 1:20-1:200

Validations

A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:108826(CALML3 Antibody) at dilution of 1:300



Immunofluorescent analysis of A431 cells using Catalog No:108826(CALML3 Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Rabbit IgG