

Calpain 7 antibody



Catalog Number: 108837

Product name

Polyclonal

Calpain 7 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

Calpain 7 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 92 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of Calpain 7 recombinant protein (Accession Number: NM_014296). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:108837(CAPN7 antibody) at dilution of 1:1000