

Carboxypeptidase A5 antibody



Catalog Number: 108850

Product name

Carboxypeptidase A5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Carboxypeptidase A5 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Positive IF detected in PC-3 cells. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of Carboxypeptidase A5 recombinant protein (Accession Number: XM_024446999). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108850(CPA5 antibody) at dilution of 1:500



Immunofluorescent analysis of PC-3 cells using Catalog No:108850(CPA5 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)