

# Carboxypeptidase M antibody



Catalog Number: 108853

## Product name

Carboxypeptidase M antibody

## Specificity

Human; other species not tested.

## Antibody description

Carboxypeptidase M Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 54 kDa

## Preparation

This antibody was obtained by immunization of Carboxypeptidase M recombinant protein (Accession Number: NM\_198320). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108853(CPM antibody) at dilution of 1:600