

CARD9 antibody



Catalog Number: 108858

Product name

CARD9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CARD9 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HEK-293 cells, HL-60 cells, human placenta tissue. Observed molecular weight by Western-blot: 50-59kd

Preparation

This antibody was obtained by immunization of CARD9 recombinant protein (Accession Number: NM_052814). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108858(CARD9



antibody) at dilution of 1:500

HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:108858(CARD9 antibody) at dilution of 1:300