

CASK antibody



Catalog Number: 108869

Product name

CASK antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CASK Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse brain tissue, mouse liver tissue. Observed molecular weight by Western-blot: 100-110kd, 50-55kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003688). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108869(CASK antibody) at dilution of 1:300