

Caspase 2/p18 antibody



Catalog Number: 108873

Product name

Caspase 2/p18 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Caspase 2/p18 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HEK-293 cells, HeLa cells, HL-60 cells, Jurkat cells, mouse liver tissue. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of Caspase 2/p18 recombinant protein (Accession Number: NM_032982). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108873(CASP2 antibody) at dilution of 1:500