

# Caspase 2/p18 antibody



Catalog Number: 108873

## Product name

Caspase 2/p18 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Caspase 2/p18 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HEK-293 cells, HeLa cells, HL-60 cells, Jurkat cells, mouse liver tissue. Observed molecular weight by Western-blot: 51 kDa

## Preparation

This antibody was obtained by immunization of Caspase 2/p18 recombinant protein (Accession Number: NM\_032982). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108873(CASP2 antibody) at dilution of 1:500