

# Caspase 2/p18 antibody



Catalog Number: 108874

## Product name

Caspase 2/p18 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Caspase 2/p18 Rabbit Polyclonal antibody.  
Positive WB detected in HEK-293 cells, Jurkat cells, mouse kidney tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 45 kDa, 32 kDa

## Preparation

This antibody was obtained by immunization of Caspase 2/p18 recombinant protein (Accession Number: NM\_032982). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108874(Caspase 2 Antibody) at dilution of



Immunofluorescent analysis of (-20oc Ethanol) fixed HEK-293 cells using Catalog No:108874(Caspase 2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)