

Caspase 2/p18 antibody



Catalog Number: 108874

Product name

Caspase 2/p18 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Caspase 2/p18 Rabbit Polyclonal antibody.
Positive WB detected in HEK-293 cells, Jurkat cells, mouse kidney tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 45 kDa, 32 kDa

Preparation

This antibody was obtained by immunization of Caspase 2/p18 recombinant protein (Accession Number: NM_032982). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108874(Caspase 2 Antibody) at dilution of



Immunofluorescent analysis of (-20oc Ethanol) fixed HEK-293 cells using Catalog No:108874(Caspase 2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)