

Caspase 4/p20/p10 antibody



Catalog Number: 108877

Product name

Caspase 4/p20/p10 antibody

Specificity

Human; other species not tested.

Antibody description

Caspase 4/p20/p10 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in SGC-7901 cells, human testis tissue. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of Caspase 4/p20/p10 recombinant protein (Accession Number: NM_001225). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:108877(CASP4 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:108877(CASP4 antibody) at dilution of 1:50 (under 25x lens)