

Caspase 6/p18/p11 antibody



Catalog Number: 108878

Product name

Caspase 6/p18/p11 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Caspase 6/p18/p11 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human lung cancer tissue, human skin tissue. Positive WB detected in Jurkat cells, HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 33-35 kDa

Preparation

This antibody was obtained by immunization of Caspase 6/p18/p11 recombinant protein (Accession Number: NM_001226). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108878(CASP6 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:108878(CASP6 antibody) at dilution of 1:50 (under 10x lens)