

# CASPR antibody



Catalog Number: 108882

## Product name

CASPR antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CASPR Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IHC detected in human cerebellum tissue. Observed molecular weight by Western-blot: 180-190 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003632). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108882(CASPR antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human cerebellum using Catalog No:108882(CASPR antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human cerebellum using Catalog No:108882(CASPR antibody) at dilution of 1:50 (under 40x lens)