

Cathepsin G antibody



Catalog Number: 108891

Product name

Cathepsin G antibody

Immunogen

[Human CTSG Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

Cathepsin G Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 29 kDa, 35 kDa

Preparation

This antibody was obtained by immunization of Cathepsin G recombinant protein (Accession Number: NM_001911). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108891(CTSG Antibody) at dilution of 1:600