

# Cathepsin G antibody



Catalog Number: 108891

## Product name

Cathepsin G antibody

## Immunogen

[Human CTSG Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

Cathepsin G Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 29 kDa, 35 kDa

## Preparation

This antibody was obtained by immunization of Cathepsin G recombinant protein (Accession Number: NM\_001911). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108891(CTSG Antibody) at dilution of 1:600