

# CBL antibody



Catalog Number: 108898

## Product name

CBL antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CBL Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells, Jurkat cells, mouse spleen tissue, Raji cells. Observed molecular weight by Western-blot: 100 kDa

## Preparation

This antibody was obtained by immunization of CBL recombinant protein (Accession Number: NM\_005188). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:108898(CBL Antibody) at dilution of 1:1000



IP Result of anti-CBL (IP:Catalog No:108898, 4ug; Detection:Catalog No:108898 1:1000) with K-562 cells lysate 3200ug.