

Calponin antibody



Catalog Number: 108906

Product name

Calponin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Calponin Rabbit Polyclonal antibody. Positive WB detected in mouse stomach tissue. Positive IHC detected in human heart tissue, human breast hyperplasia tissue, human skeletal muscle tissue. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of Calponin recombinant protein (Accession Number: NM_001299). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



mouse stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:108906(Calponin Antibody) at dilution of 1:3000



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:108906(Calponin Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog

Calponin antibody



Catalog Number: 108906

No:108906(Calponin Antibody) at dilution of 1:50 (under 40x lens)