

CAMK1D antibody



Catalog Number: 108912

Product name

CAMK1D antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CAMK1D Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 40kd

Preparation

This antibody was obtained by immunization of CAMK1D recombinant protein (Accession Number: NM_153498). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:108912(CAMK1D antibody) at dilution of 1:500