

CAMK2N1 antibody



Catalog Number: 108914

Product name

CAMK2N1 antibody

Specificity

Human; other species not tested.

Antibody description

CAMK2N1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 10 kDa

Preparation

This antibody was obtained by immunization of CAMK2N1 recombinant protein (Accession Number: BC020630). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108914(CAMK2N1 Antibody) at dilution of 1:600