

CAMTA2 antibody



Catalog Number: 108929

Product name

CAMTA2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CAMTA2 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of CAMTA2 recombinant protein (Accession Number: BC032381). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108929(CAMTA2 Antibody) at dilution of 1:700