

CANT1 antibody

Catalog Number: 108931

Product name

CANT1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CANT1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 34kd,40kd

Preparation

This antibody was obtained by immunization of CANT1 recombinant protein (Accession Number: XM_011524294). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:108931(CANT1 antibody) at dilution of 1:300