

CBX6 antibody

Catalog Number: 108948

Product name

CBX6 antibody

Specificity

Human; other species not tested.

Antibody description

CBX6 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IHC detected in human skin tissue. Observed molecular weight by Western-blot: 55-60 kDa

Preparation

This antibody was obtained by immunization of CBX6 recombinant protein (Accession Number: NM_014292). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108948(CBX6 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human skin slide using Catalog No:108948(CBX6 Antibody) at dilution of 1:50