

CCAR2 antibody



Catalog Number: 108954

Product name

CCAR2 antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:500-1:5000


IHC: 1:20-1:200

Antibody description

CCAR2 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, HeLa cells, MCF7 cells, rat brain tissue. Positive IHC detected in human breast cancer tissue, human colon cancer tissue, human pancreas cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 130 kDa

IF: 1:20-1:200

Validations




Preparation

This antibody was obtained by immunization of CCAR2 recombinant protein (Accession Number: BC018269). Purification method: Antigen Affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108954(DBC1 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:108954(DBC1 Antibody) at dilution of 1:400


Applications

ELISA, WB, IF, IHC, IP


Dilutions


CCAR2 antibody

Catalog Number: 108954



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:108954(DB1 Antibody) at dilution of 1:400



Immunofluorescent analysis of HepG2 cells, using KIAA1967 antibody Catalog No:108954 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

IP Result of anti-DB1 (IP:Catalog No:108954, 4ug; Detection:Catalog No:108954 1:1000) with mouse brain tissue lysate 4000ug.