

CCDC42 antibody



Catalog Number: 108975

Product name

CCDC42 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CCDC42 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, human brain tissue. Positive IHC detected in human brain. Observed molecular weight by Western-blot: 38kd

Preparation

This antibody was obtained by immunization of CCDC42 recombinant protein (Accession Number: NM_001158261). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108975(CCDC42 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:108975(CCDC42 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:108975(CCDC42 antibody) at dilution of 1:100 (under 40x lens)