

# CCDC42 antibody



Catalog Number: 108975

## Product name

CCDC42 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

CCDC42 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, human brain tissue. Positive IHC detected in human brain. Observed molecular weight by Western-blot: 38kd

## Preparation

This antibody was obtained by immunization of CCDC42 recombinant protein (Accession Number: NM\_001158261). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108975(CCDC42 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:108975(CCDC42 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:108975(CCDC42 antibody) at dilution of 1:100 (under 40x lens)