

CCDC56 antibody



Catalog Number: 108981

Product name

CCDC56 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CCDC56 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in Jurkat cells, human heart tissue. Observed molecular weight by Western-blot: 12kd

Preparation

This antibody was obtained by immunization of CCDC56 recombinant protein (Accession Number: NM_001040431). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108981(CCDC56 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using CCDC56 antibody Catalog No:108981 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).