

CCDC6 antibody



Catalog Number: 108982

Product name

CCDC6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CCDC6 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of CCDC6 recombinant protein (Accession Number: NM_005436). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108982(CCDC6 antibody) at dilution of 1:500



Immunofluorescent analysis of MCF-7 cells using Catalog No:108982(CCDC6 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)